

VENTILATION TUBES



MICROGEL™

Surface-Enhanced Silicone Vent Tubes






Reduces bacterial adherence by more than 90%

Myringotomy tube failure and premature extrusion can often be linked to the adherence and proliferation of bacteria on the surface of the material. Silicone substrates are especially prone to bacterial attachment because of their inherent surface morphology.¹

Microgel™ Surface-Enhanced Silicone Vent Tubes have a modified polyvinylpyrrolidone (PVP) hydrogel coating. The patented Microgel process increases the surface energy of silicone, reducing the attachment of microbial pathogens.² Bacterial adherence assays have shown that the Microgel surface enhancement reduces the occurrence of bacterial adherence by more than 90%.

Microgel™ vent tubes are available in a variety of designs with the patented, modified polyvinylpyrrolidone (PVP) hydrogel coating. This surface enhancement reduces bacterial adherence by more than 90%.



Staphylococcus epidermidis


Composition

Microgel™ is covalently bonded to the surface of silicone for uniform and long-lasting surface characteristics.

¹Karlan, M.S. Myringotomy Tube Materials: Bacterial Adhesion and Infection. *Otolaryngol Head Neck Surgery*. 1980; 88:783-795.
²Dunkirk, S.G., et al. Photochemical Coatings for the Prevention of Bacterial Colonization. *Journal of Biomaterials Applications*. 1991; 91:131-156.

Armstrong Beveled Grommets


Designed by Beverly W. Armstrong, MD, Charlotte, NC



Material	ID	IFD	Product#
Fluoroplastic, Green, 5/box	1.14	3.5	10-10030
Fluoroplastic, Green, 50/box	1.14	3.5	10-10031
Fluoroplastic, Green w/wire, 5/box	1.14	3.5	10-10035
Activent® Fluoroplastic, 5/box	1.14	3.5	10-30030
Microgel™, w/Tab, Green, 5/box	1.14	3.5	2550-1
Microgel™, w/Tab, Green, 5 pr/box	1.14	3.5	2550-2
Activent® Silicone, 5/box	1.14	3.5	10-26055
Activent® Silicone, 5 pr/box	1.14	3.5	10-46055
Firm Silicone, w/Tab, 5/box	1.14	3.5	2450-1
Firm Silicone, w/Tab, 5 pr/box	1.14	3.5	2450-2
Silicone, 5/box	1.14	3.5	10-10055
Silicone, 50/box	1.14	3.5	10-10056
H/C-Flex, 5/box	1.14	3.5	10-56032
H/C-Flex, 5 pr/box	1.14	3.5	10-66032
H/C-Flex, 50/box	1.14	3.5	10-76032

Modified Armstrong Beveled Grommets

Designed by Beverly W. Armstrong, MD, Charlotte, NC



Material	ID	IFD	Product#
Fluoroplastic, White, 5/box	1.14	3.5	1442-1
Fluoroplastic, White, 5/box	1.14	3.5	10-14242
Fluoroplastic, White, 50/box	1.14	3.5	10-14243
Fluoroplastic, White, 5pr/box	1.14	3.5	10-40024
Fluoroplastic, Green, 5/box	1.14	3.5	10-10170
Fluoroplastic, Green, 50/box	1.14	3.5	10-10171
Fluoroplastic, Green, 5 pr/box	1.14	3.5	10-40045
Activent Fluoroplastic, 5/box	1.14	3.5	10-30170
Microgel™, Blue, 5/box	1.14	3.5	10-81101
Microgel™, Blue, 5pr/box	1.14	3.5	10-81102
Microgel™, Blue, 50/box	1.14	3.5	10-81103

Armstrong V Grommets

Designed by Beverly W. Armstrong, MD, Charlotte, NC



Material	ID	IFD	Product#
H/C-Flex, 5/box	1.14	2.1	10-56033
H/C-Flex, 5 pr/box	1.14	2.1	10-66033

Baxter Beveled Buttons


Designed by David Baxter, RN, Boston, MA



Material	ID	IFD	L	Product#
Fluoroplastic, White, 5/box	1.27	3.0	2.1	1461-1
Fluoroplastic, White, 5 pr/box	1.27	3.0	2.1	1461-2
Fluoroplastic, White, 5/box	0.97	2.3	1.6	1462-1
Fluoroplastic, White, 5 pr/box	0.97	2.3	1.6	1462-2
Activent® Fluoroplastic, 5/box	1.27	3.0	2.1	3461-1
Activent® Silicone, 5/box	1.27	2.6	2.1	2663-1
Firm Silicone, 5/box	1.27	2.6	2.1	2463-1
Firm Silicone, 5pr/box	1.27	2.6	2.1	2463-2

Goode T-Grommets

Endorsed by Richard L. Goode, MD, Palo Alto, CA




Material	ID	FL	IFDs	Product#
Silicone, 5/box	1.27	8.5	2.5	10-16020
Silicone, w/ Membrane, 5/box	1.27	8.5	2.5	10-16025
Activent® Silicone, 5/box	1.27	8.5	2.5	10-26011

Hubbard Airplane Tube

Endorsed by Ron Hubbard, MD, Memphis, TN

Acknowledgement to Joseph B. Touma, MD, Huntington, WV



Material	ID	IFD	Product#
Fluoroplastic, White, 5/box	0.76	2.0	10-56170

Lindeman-Silverstein Arrow

Designed by Roger C. Lindeman, MD, Seattle, WA; and Herbert Silverstein, MD, Sarasota, FL



Material	ID	IFD	L	Product#
Fluoroplastic, Green, 5/box	1.14	3.0	4.7	10-21000

Moretz Tab Tubes


Endorsed by William H. Moretz, Jr., MD, Augusta, GA



Material	ID	IFDs	Product#
Microblue® Titanium, w/Tab, 5/box	1.27	1.0	1464-1
Microblue® Titanium, w/Tab, 5 pr/box	1.27	1.0	1464-2
Tytan® Titanium, 5/box	1.02	1.0	10-56104
Tytan® Titanium, 5/box	1.27	1.0	10-56107
Tytan® Titanium, 5 pr/box	1.27	1.0	10-66107
Tytan® Titanium, 5/box	0.76	1.0	10-56108
Tytan® Titanium 5 pr/box	0.76	1.0	10-66108
Tytan® Titanium, w/ Wire, 5/box	1.27	1.0	10-56113
Tytan® Titanium, w/ Wire, 5 pr/box	1.27	1.0	10-66113
Activent® Fluoroplastic, 5/box	1.27	1.0	10-36177
Fluoroplastic, white, 5/box	1.27	1.0	10-56177
Fluoroplastic, white, 5 pr/box	1.27	1.0	10-66177
H/C-Flex, 5/box	1.27	0.75	10-56705

Reuter Bobbins

Designed by S. Harold Reuter, MD, Houston, TX



Material	ID	IFD	IFDs	Product#
Fluoroplastic, White, w/o Holes, 5/box	1.02	2.7	1.1	1412-1
Fluoroplastic, White, w/o Holes, 5 pr/box	1.02	2.7	1.1	1412-2
Fluoroplastic, White, w/o Holes, 5/box	1.14	2.7	1.1	1413-1
Fluoroplastic, White, w/o Holes, 5 pr/box	1.14	2.7	1.1	1413-2
Fluoroplastic, White, w/ Holes, 5/box	1.14	2.7	1.1	1425-1
Fluoroplastic, White, w/ Holes, 5 pr/box	1.14	2.7	1.1	1425-2
Fluoroplastic, Blue, w/ Holes, 5/box	1.14	2.6	1.0	10-10101
Fluoroplastic, Blue, w/ Holes, 50/box	1.14	2.6	1.0	10-10103
Fluoroplastic, Blue, w/o Holes, 5/box	1.14	2.6	1.0	10-10201
Fluoroplastic, Blue, w/o Holes, 50/box	1.14	2.6	1.0	10-10203
Fluoroplastic, Blue, w/o Holes, 5 pr/box	1.14	2.6	1.0	10-10202
Fluoroplastic, Blue, w/ Holes & Wire, 5/box	1.14	2.6	1.0	10-10301
Fluoroplastic, White, w/o Holes, 5/box	1.02	2.6	1.0	10-56162
Activent® Fluoroplastic w/o Holes, 5/box	1.02	2.6	1.0	10-36162
Activent® Fluoroplastic, w/ Holes, 5/box	1.14	2.6	1.0	10-37000
Activent® FLPL w/o Holes, 5/box	1.14	2.6	1.0	10-37005
Stainless Steel, Matte Finish, w/ Holes, 5/box	1.02	2.7	1.0	14181
Stainless Steel, w/ Holes 5/box	1.27	2.7	1.0	14051
Stainless Steel, w/ Holes & Wire, 5/box	1.14	2.6	1.0	10-27040
Stainless Steel, w/ Holes, 5/box	1.02	2.7	1.0	1440-1
Stainless Steel, w/ Holes, 5 pr/box	1.02	2.7	1.0	1440-2
Stainless Steel, w/ Holes, 5/box	1.27	2.6	0.5	10-27010
Stainless Steel, w/ Holes, 5/box	1.14	2.6	1.0	10-27015
Stainless Steel, w/ Holes & Wire, 5 pr/box	1.27	2.7	1.0	1406-2
Stainless Steel, w/ Holes & Wire, 50/box	1.14	2.6	1.0	10-27016
Stainless Steel, w/ Holes, 5 pr/box	1.14	2.6	1.0	10-40041
Titanium, w/o Holes, 5/box	1.27	2.6	1.0	10-27070
Titanium, w/ Holes, 5/box	1.02	2.7	1.0	1611-1
Titanium, w/ Holes, 5/box	1.02	2.7	1.0	1611-2
Titanium, w/ Holes, 5/box	1.14	2.6	1.0	10-27060
Titanium, w/ Holes, 5/box	1.27	2.6	0.5	10-27065
Tytan® Titanium, w/ Holes, 5/box	1.02	2.7	1.0	10-56110
Tytan® Titanium, w/ Holes, 5 pr/box	1.02	2.7	1.0	10-66110
Tytan® Titanium, w/o Holes, 5/box	1.14	2.7	1.0	10-56112
Tytan® Titanium, w/o Holes, 5 pr/box	1.14	2.7	1.0	10-66112
Microblue, Tytan® Titanium w/o Holes, 5/box	1.27	2.7	1.0	1622-1 or 2
Stainless Steel, Microblue, Matte Finish, Titanium, w/ Holes, 5/box	1.14	2.6	1.0	10-27017
Stainless Steel, Microblue, Matte Finish, Titanium, w/ Holes, 50/box	1.14	2.6	1.0	10-27018
Stainless Steel, Microblue, Matte Finish, Titanium, w/ Holes, 5 pr/box	1.14	2.6	1.0	10-40004

Baxter Beveled T-Grommets

Designed by David Baxter, RN,
Boston, MA



Material	ID	FL	L	Product#
Firm Silicone, 5/box	1.27	9.5	3.0	2465-1
Firm Silicone, 5 pr/box	1.27	9.5	3.0	2465-2

Beveled Bobbins



Material	ID	Product#
Fluoroplastic, White, 5/box	1.14	10-56179
Fluoroplastic, White, 5pr/box	1.14	10-66179


Cohen T-Grommets

Designed by Richard Cohen, MD,
Grants Pass, OR



Material	ID	FL	Product#
Firm Silicone, 5/box	1.14	9.1	2470-1
Firm Silicone, 5 pr/box	1.14	9.1	2470-2
Microgel™, Silicone, Green, 5/box	1.14	9.1	2569-1
Microgel™, Silicone, Green, 5 pr/box	1.14	9.1	2569-2

Donaldson Type Tubes



Material	ID	IFD	IFDis	Product#
Fluoroplastic, Blue, 5/box	1.14	2.3	0.9	10-15101
Fluoroplastic, Blue, 5pr/box	1.14	2.3	0.9	10-15102
Fluoroplastic w/ Wire, Blue, 5/box	1.14	2.3	0.9	10-15201
Activent Fluoroplastic, 5/box	1.14	2.3	0.9	10-33000
Microgel™, w/ Tab, Blue, 5/box	1.14	2.3	0.9	10-85501
Microgel, w/ Tab, Blue, 5 pr/box	1.14	2.3	0.9	10-85502
Activent® Silicone, 5/box	1.14	2.3	0.9	10-26020
Activent® Silicone, 5pr/box	1.14	2.3	0.9	10-46020
Firm Silicone, 5/box	1.14	2.3	0.9	2415-1
Firm Silicone, 5 pr/box	1.14	2.3	0.9	2415-2
Firm Silicone, w/ Tab, 5/box	1.14	2.3	0.9	2416-1
Firm Silicone, w/ Tab, 5 pr/box	1.14	2.3	0.9	2416-2
Silicone, w/Membrane, 5/box	1.14	2.3	0.9	10-12015
Silicone, w/ Tab, 5/box	1.14	2.3	0.9	10-13015
Silicone, 5/box	1.14	2.3	0.9	10-13020
Silicone, 50/box	1.14	2.3	0.9	10-13021
Silicone, 5 pr/box	1.14	2.3	0.9	10-40009
C-Flex, 5/box	1.14	2.3	0.76	10-56010
C-Flex, 5 pr/box	1.14	2.3	0.76	10-66010

Paparella Type Vent Tubes


Designed by Michael M. Paparella, MD,
Minneapolis, MN



Material	ID	IFD	IFDis	Product#
Microgel™ w/ Notch & Tab, Green, 5/box	1.02	2.5	1.0	2544-1
Microgel™ w/ Notch & Tab, Green, 5 pr/box	1.02	2.5	1.0	2544-2
Microgel™ w/ Notch & Tab, Green, 50/box	1.02	2.5	1.0	2544-3
Microgel™ w/ Notch & Tab, Blue, 5/box	1.14	2.1	1.1	10-82201
Microgel™ w/ Notch & Tab, Blue, 5pr/box	1.14	2.1	1.1	10-82202
Microgel™ w/ Notch & Tab, Blue, 50/box	1.14	2.1	1.1	10-82203
Activent® Silicone w/ Notch & Tab, 5/box	1.14	2.4	1.1	10-26001
Activent® Silicone w/ Notch & Tab, 5 pr/box	1.14	2.4	1.1	10-46001
Activent® Silicone w/ Notch & Tab, 5/box	1.52	4.45	1.15	10-26045
Firm Silicone w/ Notch & Tab, 5/box	1.02	2.5	1.0	2444-1
Firm Silicone w/ Notch & Tab, 5 pr/box	1.02	2.5	1.0	2444-2
Firm Silicone w/ Notch, 5/box	1.02	2.5	1.0	2445-1
Firm Silicone w/ Notch, 5 pr/box	1.02	2.5	1.0	2445-2
Firm Silicone w/ Notch & Tab, 5/box	1.27	4.3	0.9	2446-1
Firm Silicone w/ Notch & Tab, 5 pr/box	1.27	4.3	0.9	2446-2
Biolite-coated Silicone, 5/box	1.14	2.4	1.1	10-11015
Silicone w/ Membrane, 5/box	1.14	4.3	.85	10-12010
Silicone w/ Notch & Tab, 5/box	1.14	2.4	1.1	10-25001
Silicone, w/ Notch & Tab, 50/box	1.14	2.4	1.1	10-25002
Silicone, 5 pr/box	1.14	2.4	1.1	10-40006
Silicone, 50 pr/box	1.14	2.4	1.1	10-40019
Silicone w/ Notch, 5/box	1.14	2.4	1.1	10-25005
Silicone w/ Notch, 5/box	1.52	4.45	1.15	10-25040
Silicone w/ Notch & Tab, 5/box	1.52	4.45	1.15	10-25045
Silicone w/ Notch & Tab, 50/box	1.52	4.45	1.15	10-25046
C-Flex w/ Tab, 5/box	1.02	2.5	1.0	10-56005
C-Flex, w/ Tab, 5 pr/box	1.02	2.5	1.0	10-66005

Pediatric Ventilation Tubes

10-56101, 10-56102, 10-56159 Endorsed by
M. Coyle Shea, MD, Memphis, TN
10-66108 Endorsed by
William H. Moretz, Jr., MD, Augusta, GA
1431 Designed by N. Shah, FRCS, DLO, London, England



Material	ID	IFD	Product#
Microblue Titanium, 5/box	0.76	2.0	1619-1
Tiny Tytan® Titanium, 5/box	0.76	1.5	10-56101
Tytan® Titanium, 5/box	0.76	2.0	10-56102
Tiny Tytan® Titanium, Moretz tab, 5/box	0.76	1.5	10-56108
Tiny Tytan® Titanium, Moretz tab, 5 pr/box	0.76	1.5	10-66108
Fluoroplastic, White, Mini Shah, 5/box	0.76	1.5	1431-1
Fluoroplastic, White, 5/box	0.76	2.0	10-56159
Fluoroplastic, White, 5 pr/box	0.76	2.0	10-66159
Activent® Fluoroplastic, 5/box	0.76	2.0	10-36159

Pope Beveled Grommets


Designed by Thad Pope, MD,
Crawfordsville, IN



Material	ID	IFD	IFDis	Product#
Fluoroplastic, Green, 5/box	1.14	3.6	1.3	10-25004
Fluoroplastic, Green, 5 pr/box	1.14	3.6	1.3	10-40014
Fluoroplastic, Green, 50/box	1.14	3.6	1.3	10-25007
Fluoroplastic w/ Wire Green, 5/box	1.14	3.6	1.3	10-25006
Activent® Fluoroplastic, 5/box	1.14	3.6	1.3	10-35004
H/C-Flex, 5/box	1.14	3.6	1.0	10-56036
H/C-Flex, 5 pr/box	1.14	3.6	1.0	10-66036
H/C-Flex, w/ Tail, 5/box	1.14	3.6	1.0	10-56039
H/C-Flex, w/ Tail, 5 pr/box	1.14	3.6	1.0	10-66039

Rube T-Grommets

Designed by William Rubenstein, MD,
Houston, TX



Material	ID	FL	Product#
Firm Silicone, 5/box	1.14	6.5	2462-1

Shah Type Ventilation Tubes

Designed by N. Shah, FRCS, DLO,
London, England



Material	ID	IFDis	Product#
Fluoroplastic, Green, 5/box	1.14	1.58	10-28030
Fluoroplastic, Green, w/ Wire, 5/box	1.14	1.58	10-28035
Fluoroplastic, Green, w/ Wire, 5/box	1.14	2.2	10-56174
Activent® Fluoroplastic, 5/box	1.14	1.58	10-38030
Fluoroplastic, White, w/ Wire, 5 pr/box	1.14	1.75	1429-2
Fluoroplastic, White, Mini, 5/box	0.76	1.00	1431-1

Sheehy Type Collar Buttons

Designed by James L. Sheehy, MD,
Los Angeles, CA



Material	ID	IFD	IFDis	Product#
Fluoroplastic, White, 5/box	1.27	3.0	1.52	1414-1
Fluoroplastic, White, 5/box	1.27	3.0	1.52	1414-2
Fluoroplastic, White, w/ Wire, 5/box	1.27	3.0	1.52	1416-1
Fluoroplastic, White, 5/box	1.27	3.1	1.5	10-56160
Fluoroplastic, White, 5 pr/box	1.27	3.1	1.5	10-66160
Fluoroplastic, Blue, 5/box	1.27	3.0	1.55	10-13301
Fluoroplastic, Blue, 50/box	1.27	3.0	1.55	10-13303
Fluoroplastic, Blue, 5 pr/box	1.27	3.0	1.55	10-13302
Fluoroplastic, Blue w/ Wire, 5/box	1.27	3.0	1.55	10-13201
Activent® Fluoroplastic, 5/box	1.27	3.0	1.55	10-38040
Activent® Silicone, 5/box	1.27	2.9	1.3	10-26145
Activent® Silicone, 50/box	1.27	2.9	1.3	10-26146
Activent® Silicone, 5 pr/box	1.27	2.9	1.3	10-46145
Silicone, w/ Membrane, 5/box	1.27	2.8	1.32	10-12050
Silicone, 5/box	1.27	2.9	1.3	10-28145
Silicone, 50/box	1.27	2.9	1.3	10-28146
H/C-Flex, 5/box	1.27	3.1	1.5	10-56037
H/C-Flex, 5 pr/box	1.27	3.1	1.5	10-66037
Tytan®, Titanium, 5/box	1.27	3.1	1.5	10-56114
Microgel™, Blue, 5/box	1.27	2.9	1.3	10-83301
Microgel™, Blue, 5pr/box	1.27	2.9	1.3	10-83302
Microgel™, Blue, 50/box	1.27	2.9	1.3	10-83303

Shepard Grommets



Material	ID	IFD	Product#
Fluoroplastic, Blue, w/ Wire, 5/box	1.02	2.4	10-16501
Fluoroplastic, Blue, w/ Wire, 5 pr/box	1.02	2.4	10-16502
Fluoroplastic, Blue, w/ Wire, 50/box	1.02	2.4	10-16503
Fluoroplastic, Blue, 5/box	1.02	2.4	10-16301
Microgel™, w/Tab, Blue, 5/box	1.02	2.3	10-86601
Microgel™, w/Tab, Blue, 5 pr/box	1.02	2.3	10-86602
Fluoroplastic, Green, w/ Tail, 5/box	1.14	2.4	10-28010
Fluoroplastic, Green, w/ Tail, 5 pr/box	1.14	2.4	10-40055
Fluoroplastic, Green, w/ Tail, 50/box	1.14	2.4	10-28011
Fluoroplastic, Blue, w/ Wire, 5/box	1.14	2.4	10-16201
Fluoroplastic, Blue, w/ Wire, 50/box	1.14	2.4	10-16203
Fluoroplastic, w/ Tab, 5/box	1.14	2.4	10-28015
Firm Silicone, w/ Tab, 5/box	1.02	2.3	2420-1
Firm Silicone, w/ Tab, 5 pr/box	1.02	2.3	2420-2
Fluoroplastic, Blue, w/ Tab, 5/box	1.14	2.4	10-16101
Fluoroplastic, Blue, w/ Tab, 5 pr/box	1.14	2.4	10-16102
Fluoroplastic, Blue, w/ Tab, 50/box	1.14	2.4	10-16103
Fluoroplastic, White, 5/box	1.14	2.5	10-56154
Fluoroplastic, White, 5 pr/box	1.14	2.5	10-66154
Fluoroplastic, White, 5/box	1.02	2.5	10-56155
Fluoroplastic, White, 5 pr/box	1.02	2.5	10-66155
Fluoroplastic, White, w/ Wire, 5/box	1.14	2.5	10-56157
Activent® Fluoroplastic, w/ Wire, 5/box	1.14	2.4	10-38012
Activent®, Silicone, w/ Tail, 5/box	1.14	2.3	10-26025
Activent®, Silicone, w/ Tail, 5 pr/box	1.14	2.3	10-46025
Silicone, w/ Tail, 5/box	1.14	2.3	10-28025
Silicone, w/ Tail, 50/box	1.14	2.3	10-28026

Soileau Tytan® Titanium

Endorsed by James S. Soileau, MD, Gonzales, LA



Material	ID	IFD	Product#
Tytan® Titanium, 5/box	1.27	2.75	10-56111
Tytan® Titanium, 5 pr/box	1.27	2.75	10-66111
Tytan® Titanium, 50/box	1.27	2.75	10-76111

Spoon Bobbins



Material	ID	IFD	Product#
Tytan® Titanium, 5 pr/box	1.27	2.7	10-66106
Tytan® Titanium, 5/box	1.02	2.7	10-56109
Tytan® Titanium, 5/box	1.27	2.7	10-56106
Stainless Steel, w/ Holes, 5 pr/box	1.02	2.7	1400-2

Armstrong Beveled Straight Shanks

Designed by Beverly W. Armstrong, MD, Charlotte, NC



Material	ID	IFD	L	Product#
Fluoroplastic, Green, 5/box	1.14	3.6	10.0	10-10150
Fluoroplastic, Green, 5/box	1.14	3.6	8.0	10-10005
Fluoroplastic, Green, 5/box	1.14	3.6	7.0	10-10000
Fluoroplastic, Green, 5 pr/box	1.14	3.6	7.0	10-40022
Silicone, 5/box	1.14	3.6	12.0	10-10050
Activent® Silicone, 5/box	1.14	3.6	7.0	10-26140
Silicone, 5/box	1.14	3.6	7.0	10-10040
Silicone, 5/box	1.14	3.6	10.0	10-10045

Armstrong V Tubes

Designed by Beverly W. Armstrong, MD, Charlotte, NC



Material	ID	L	Product#
Fluoroplastic, White, 5/box	1.14	7.5	10-56120
Fluoroplastic, White, 5 pr/box	1.14	7.5	10-661202
Fluoroplastic, White, 5/box	1.14	4.5	10-56121
Fluoroplastic, White, 5 pr/box	1.14	4.5	10-661212
Activent® Fluoroplastic, 5/box	1.14	4.5	10-36121
H/C-Flex, 5/box	1.14	7.5	10-56030
H/C-Flex, 5/box	1.14	4.5	10-56031

Goode T-Tubes®

Designed by Richard L. Goode, MD, Palo Alto, CA



Material	ID	FL	L	Product#
Microgel™, Silicone, Blue, 5/box	1.14	9.8	12.0	10-87701
Microgel™, Silicone, Blue, 5 pr/box	1.14	9.8	12.0	10-87702
Activent® Silicone, 5/box	1.14	9.5	12.0	10-26010
Activent® Silicone, Blue, 50/box	1.14	9.5	12.0	10-37703
Biolite-coated Silicone, 5/box	1.14	9.5	12.0	10-11030
Silicone, 5/box	1.14	9.5	12.0	10-16010
Silicone, 50/box	1.14	9.5	12.0	10-16011
Silicone, 5 pr/box	1.14	9.5	12.0	10-40001
Firm Silicone, 5/box	1.14	9.5	12.0	10-16035
Firm Silicone, 5 pr/box	1.14	9.5	12.0	2471-2
Firm Silicone, 5/box	1.14	6.0	6.0	2475-1
Firm Silicone, 5 pr/box	1.14	6.0	6.0	2475-2
C-Flex, 5/box	1.14	9.0	12.0	10-56015
C-Flex, 5/box	1.14	6.5	6.0	10-56016
C-Flex, 5 pr/box	1.14	6.5	6.0	10-66016
C-Flex, 5/box	1.14	8.0	9.0	10-56019

Hoffman Permanent Vent Tubes



Material	ID	L	Product#
Titanium, w/HA & Fluoroplastic, w/ Liner, 1/box	1.65	8.3	10-28080
Fluoroplastic, Inner Liner Only, 5/box	1.25	13.4	10-28081

Permanent Aeration Tubes



Material	ID	FL	L	Product#
Silicone, Fluoroplastic Stylet, 5/box	1.14	3.5	12.0	10-28065
Silicone, 5/box	1.14	3.5	7.0	10-28070

Semi-Permeable Membrane Tubes

Silicone w/ Fluoroplastic Membrane



Material	ID	IFD	IFDis	Product#
Donaldson, 5/box	1.14	2.3	0.9	10-12015
Goode T-Grommet, Silicone, 5/box	1.27	8.5	2.5	10-16025
Paparella, 5/box	1.14	4.3	0.9	10-12010
Sheehy, 5/box	1.27	2.8	1.32	10-12050

Kartush Tympanic Membrane (TM) Patcher

Designed by Jack M. Kartush, MD, Farmington Hills, MI



Material	DIA	CD	Product#
Silicone, 1/box	3.0	0.9	14-16000
Silicone, 1/box	5.0	1.1	14-16010
Silicone, 1/box	7.0	1.3	14-16020
Silicone	10.0	1.4	14-16030

Touma Beveled T-Grommet


Endorsed by Joseph B. Touma, MD,
Huntington, WV



Material	ID	IFD	Product#
Activent® Silicone, w/ Tab, 5/box	1.27	5.0	10-26014

Touma T-Type Grommets

Endorsed by Joseph B. Touma, MD,
Huntington, WV



Material	ID	IFD	Product#
Activent® Silicone, 5/box	1.14	5.0	10-26012
C-Flex, 5/box	1.14	5.0	10-56012
C-Flex, 5 pr/box	1.14	5.0	10-66012

Umbrella Ventilation Tubes



Material	ID	L	Product#
Silicone, w/ Tab, 5/box	1.02	2.2	10-16030
Activent® Silicone, w/ Tab, 5/box	1.02	2.2	10-26030
Microgel™, w/ Tab, Blue, 5/box	1.02	2.2	10-88801
Microgel™, w/ Tab, Blue, 5pr/box	1.02	2.2	10-88802

Venturi™ Style Ventilation Tubes



Material	ID	IFD	IFDis	Product#
Titanium, Reuter Bobbin, 5/box	1.02	2.5	1.0	10-70011
Titanium, Reuter Bobbin, 5/box	1.27	2.5	1.0	10-70031
Titanium, Reuter Bobbin, 5 pr/box	1.27	2.5	1.0	10-70032
Microblue, Titan Reuter Bobbin, 5/box	1.27	2.5	1.0	10-70051
Microblue, Titan Reuter Bobbin, 5 pr/box	1.27	2.5	1.0	10-70052

V-T Ventilation Tube



Material	ID	IFD	IFDis	Product#
C-Flex, 5/box	1.27	4.5	1.3	10-02341
C-Flex, 50/box	1.27	4.5	1.3	10-64403
Microgel™, Blue, 5/box	1.27	4.5	1.3	10-84401
Microgel™, Blue, 5pr/box	1.27	4.5	1.3	10-84402
Microgel™, Blue, 50/box	1.27	4.5	1.3	10-84403

Modified Goode T-Tubes



Material	ID	FL	L	Product#
Activent® Silicone, 5/box	1.27	9.5	4.5	10-26040
Silicone, 5 pr/box	1.27	5.8	6.7	2477-1
Silicone, 5 pr/box	1.27	5.8	6.7	2477-2
Firm Silicone, 5/box	1.27	9.5	4.5	10-16040
Firm Silicone, 50/box	1.27	9.5	4.5	10-16041
Firm Silicone, 5 pr/box	1.27	9.5	4.5	10-40002

Per-Lee Tube

Designed by John H. Per-Lee, MD,
Atlanta, GA




Material	ID	FL	L	Product#
Silicone, 60°, 5/box	1.5	7.8	15.0	10-25020

Plain Straight Shanks




Material	ID	IFD	L	Product#
Fluoroplastic, White, 5/box	1.14	2.3	12.0	1435-1
Fluoroplastic, Blue, 5/box	1.14	2.75	7.0	10-19301
Fluoroplastic, Blue, 5/box	0.89	2.75	7.0	10-19101
Fluoroplastic, White, 5/box	0.89	2.1	7.0	10-56165
Fluoroplastic, White, 5 pr/box	0.89	2.1	7.0	10-66165
Fluoroplastic, Blue, 5/box	1.14	2.2	7.0	10-19201
Fluoroplastic, White, 5/box	0.89	2.1	12.0	10-56166
Fluoroplastic, White, 5/box	1.14	2.1	7.0	10-56167
Activent® Fluoroplastic, 5/box	1.14	2.75	7.0	10-38050
Activent® Fluoroplastic, 5/box	0.89	2.1	7.0	10-36165

Siegel T-Tubes




Material	ID	FL	L	Product#
Silicone, 5/box	1.14	4.0	7.0	2481-1
Silicone, 5 pr/box	1.14	4.0	7.0	2481-2

T-Type Vent Tubes



Material	ID	FL	L	Product#
C-Flex, 5/box	1.14	4.8	6.0	10-56018
C-Flex, 5 pr/box	1.14	4.8	6.0	10-66018


Dimensions



All dimensions are approximate and subject to normal manufacturing variance.

All measurements are listed as millimeters.


ID = Inner Diameter
OFD = Outer Flange Diameter
IFD = Inner Flange Diameter
IFDis = Interflange Distance
L = Length
FL = Flange Length
DIA = Diameter



All dimensions are approximate and subject to normal manufacturing variance.

All measurements are listed as millimeters.

ID = Inner Diameter
OFD = Outer Flange Diameter
IFD = Inner Flange Diameter
IFDis = Interflange Distance
L = Length
FL = Flange Length
DIA = Diameter



All dimensions are approximate and subject to normal manufacturing variance.

All measurements are listed as millimeters.

ID = Inner Diameter
OFD = Outer Flange Diameter
IFD = Inner Flange Diameter
IFDis = Interflange Distance
L = Length
CD = Center Diameter
DIA = Diameter

For further information, please call Medtronic ENT at 800-874-5797 or 904-296-9600. You may also consult our website at www.MedtronicENT.com.

Medtronic ENT

Medtronic USA, Inc.
6743 Southpoint Drive North
Jacksonville, FL 32216
USA
www.MedtronicENT.com
Toll free: (800) 874-5797
Fax: (800) 678-3995

International Telephone Numbers

Australia: 1-800-668-670	Korea: 82-2-3404-3600
Belgium: 32-2456-09-00	Lebanon: 961-1-370-670
Canada: 1-800-217-1617	Luxembourg: 32-2456-09-09
China: 86-21-50800998	Netherlands: 31-45-566-8000
France: 33-470-679-800	Poland: 48-22-465-6942
Germany: 49-211-5293-209	Singapore: 65-6776-6255
Hong Kong: 852-2919-1312	Spain: 34-91-625-05-40
India: 91-22-26836733	UK: 44-1923-212213
Italy: 39-02-24137-324	USA: 1-904-296-9600
Japan: 81-3-6430-2011	

This literature is intended for the exclusive use of physicians. Rx only.

©2007, Medtronic, Inc.

*Registered mark of Medtronic, Inc. ™ Trademark of Medtronic, Inc

Biolite® is a registered trademark of Carbomedics, Inc.

Microblue® is a registered mark of Material Technologies, Inc.

890149 2008-027 7.07 Rev. 3



Medtronic

Alleviating Pain · Restoring Health · Extending Life